

# A Fast, Simple, and Affordable Technique to Measure Oxygen Consumption in Living Zebrafish Embryos

Julie Somkhit,<sup>1,\*</sup> Roxane Loyant,<sup>2,\*</sup> Alexandre Brenet,<sup>1</sup> Rahma Hassan-Abdi,<sup>1</sup> Constantin Yanicostas,<sup>1</sup> Mathieu Porceddu,<sup>2</sup> Annie Borgne-Sanchez,<sup>2</sup> and Nadia Soussi-Yanicostas<sup>1</sup>

## Abstract

In all animal species, oxygen consumption is a key process that is partially impaired in a large number of pathological situations and thus provides informative details on the physiopathology of the disease. In this study, we describe a simple and affordable method to precisely measure oxygen consumption in living zebrafish larvae using a spectrofluorometer and the MitoXpress Xtra Oxygen Consumption Assay. In addition, we used zebrafish larvae treated with mitochondrial respiratory chain inhibitors, antimycin A or rotenone, to verify that our method enables precise and reliable measurements of oxygen consumption.

**Keywords:** respiration, zebrafish larvae, oxygen consumption

## Introduction

RESPIRATION IS A KEY cellular function, which is altered in a great number of disorders and animal disease models, thus providing a meaningful window to interrogate the physiology of living cells and organisms.<sup>1</sup> Automated devices allowing precise measurement of both mitochondrial respiration (oxidative phosphorylation) and glycolysis have been developed in recent years, but their price and cost of usage restrain their use on a daily basis. To overcome this limitation, we describe here a simple and affordable method that allows precise measurement of oxygen consumption in living zebrafish larvae, using a simple spectrofluorometer and the MitoXpress Xtra Oxygen Consumption Assay (Agilent), a quantitative test relying on a fluorescent probe, which is quenched in a dose-dependent manner by O<sub>2</sub> molecules (Fig. 1A).

## Results


To assess the ability of this method to detect fine changes in O<sub>2</sub> consumption in living zebrafish larvae, we measured O<sub>2</sub> consumption of pools of seven 5 days post-fertilization wild-type individuals treated with either vehicle (1% dimethyl

---

<sup>1</sup>NeuroDiderot, Inserm U1141, Université de Paris, Paris, France.

<sup>2</sup>Mitologics SAS, Faculté de Médecine, Créteil, France.

\*These authors contributed equally to this work.



**FIG. 1.** Real-time measurement of oxygen consumption in 5 dpf zebrafish larvae. **(A)** Scheme of the experimental setup used to measure oxygen consumption *in vivo* in zebrafish larvae, treated or not with MRC inhibitors, with an example of a filled 96-well microplate. Pools of seven 5 dpf larvae treated with MRC inhibitors or untreated (control) were transferred to a 96-well microplate (Greiner Bio-One International) and then incubated in the presence of the MitoXpress Xtra probe, (Supplementary Data) a fluorescent molecule whose fluorescence is quenched in a dose-dependent manner by  $O_2$  through direct molecular interactions, thus making fluorescence intensity inversely proportional to the amount of free  $O_2$  dissolved in the medium. **(B)** Real-time measurement of the oxygen consumption in larvae exposed to MRC inhibitors or untreated, and blank (E3 medium without larva). **(C)** Quantification of the oxygen consumption of larvae exposed to either 2  $\mu$ M rotenone or 10  $\mu$ M AA. Whiskers represent the standard error of the mean. Rotenone ( $n=36$ ), AA ( $n=35$ ), control ( $n=18$ ), and blank ( $n=24$ ), one-way ANOVA with Tukey's Multiple Comparison Test: \*\*\* $p<0.0001$ . AA, antimycin A; ANOVA, analysis of variance; dpf, days post-fertilization; MRC, mitochondrial respiratory chain.

sulfoxide [DMSO]), 2  $\mu$ M rotenone, a specific inhibitor of the mitochondrial respiratory chain (MRC) complex I,<sup>2</sup> or 10  $\mu$ M antimycin A (AA), a specific inhibitor of the MRC complex III (Fig. 1A).<sup>3,4</sup> As previously reported, we found that 1% DMSO concentration had no effect on larval development.<sup>5</sup>

Results showed that oxygen consumption of larvae exposed to 2  $\mu$ M rotenone and 10  $\mu$ M AA was 31% and 6%, respectively, of that observed in non-treated controls (Fig. 1B, C, Supplementary Data), indicating that the method described allows reliable measurements of oxygen consumption. Thus, this simple and direct method can be used to measure respiration *in vivo* in zebrafish intoxication and disease models.

#### Disclosure Statement

No competing financial interests exist.

**Funding Information**

This study was supported by a grant from the Agence Nationale de la Recherche (ANR-16-CE18-0010 MITOXDRUGS).

**Supplementary Material**

Supplementary Data

**References**

1. Frazier AE, Thorburn DR, Compton AG. Mitochondrial energy generation disorders: genes, mechanisms, and clues to pathology. *J Biol Chem* 2019;294:5386–5395.
2. Lindahl PE, Öberg KE. The effect of rotenone on respiration and its point of attack. *Exp Cell Res* 1961;23:228–237.
3. Huang LS, Cobessi D, Tung EY, Berry EA. Binding of the respiratory chain inhibitor antimycin to the mitochondrial Bc1 complex: a new crystal structure reveals an altered intramolecular hydrogen-bonding pattern. *J Mol Biol* 2005;351:573–597.
4. Rieske JS. Inhibitors of respiration at energy-coupling site 2 of the respiratory chain. *Pharmacol Ther* 1980;11:415–450.
5. Alavi Naini SM, Yanicostas C, Hassan-Abdi R, Blondeel S, Bennis M, Weiss RJ, *et al.* Surfen and oxalyl surfen decrease tau hyperphosphorylation and mitigate neuron deficits in vivo in a vertebrate model of tauopathy. *Transl Neurodegener* 2018;7:6.

Address correspondence to:  
*Nadia Soussi-Yanicostas*  
*NeuroDiderot*  
*NeuroDiderot, INSERM U1141*  
*Hôpital Robert Debré*  
*48 Boulevard Sérurier*  
*Paris F-75019*  
*France*

*E-mail:* nadia.soussi@inserm.fr